

Thorlabs Scope Wiring


Steps to integrate a Thorlabs ECU1 or ECU2 (BScope 2) into a ScanImage 2016 system:

- Wire the Thor ECU and ScanImage 2016 system
- Start ScanImage
- Configure ScanImage
 - Integrate ThorLabs ECU and/or BScope 2 into ScanImage
 - Configure Thorlabs ECU Scanner
 - Configure Thorlabs B-Scope 2
 - Configure Resonant Scanning System

Wire the Thor ECU and ScanImage 2016 system

Use the image below as a guide. Items in black apply to both ECU1 and ECU2/BScope 2 while items in red apply only to the ECU2/BScope 2.

- Loop the ECU's line clock output terminal back to the ECU's line clock input terminal
- Connect the Thor ECU to the scanner system
 - For ECU1, connect the single DVI connector to the resonant/galvo scanner assembly
 - For ECU2/BScope 2, connect all cables for PMT and scanner controls. There will be a 25 pin D-sub for each of the two parallel scan paths and two round cables for each PMT
- If using a BScope 2, connect the MCM5000 controller. For each stage axis (X, Y, fine Z, coarse Z, rotation) there will be two D-sub cables, one for the motor and one for the encoder. There will also be three round mirror control cables
- Connect the 68pin DAQ connector to the scanner control DAQ board of the ScanImage 2016 system
- Connect the ECU's USB cable to the ScanImage 2016 PC. Windows will automatically download the drivers (requires an Internet connection) and register a USB serial port.
- If using a BScope 2, also connect the USB cable to the MCM5000 controller. An unknown device should appear in the device manager.
 - Use device manager to update the device driver and choose the option to manually select driver location.
 - Locate the BScope_ControlV3 driver package from Thorlabs.
 - After installing, a second USB serial port should appear in the device manager.



Start ScanImage

Start ScanImage 2016 and select the option to create a new machine data file.

Configure ScanImage

The ScanImage Machine Configuration Editor will be started. Beginning with the [General ScanImage Settings](#) section, enter settings information appropriately.


As sections are configured, modify them appropriately then press Next. When the Additional Components section is reached, you can press Finish or Finish and Run ScanImage to save configured information to the given machine data file.

Machine configuration can be updated at a later time by selecting the [Machine Configuration...](#) option on the File Menu on the Main Controls panel.

Integrate ThorLabs ECU and/or BScope 2 into ScanImage

To configure ThorLabs ECU and/or BScope 2, they must first be added within the ScanImage Machine Configuration Editor.

In the [General ScanImage Settings](#) section, add "Thorlabs ECU 1" and/or "Thorlabs B-Scope 2" to the Additional Components and Plugins table, depending on your hardware.




Additional Component	Delete
Thorlabs ECU 1	X
Thorlabs B-Scope 2	X
	X

The ScanImage Machine Configuration Editor will add settings sections, accordingl

A [Thorlabs ECU Scanner Settings](#) section is added if "Thorlabs ECU 1" is selected as an additional component; and a [ThorLabs BScope2 Setup](#) section is added if "Thorlabs B-Scope 2" is selected as an additional component.

Configure Thorlabs ECU Scanner

Select the ThorLabs ECU Scanner section In the ScanImage Machine Configuration Editor.



ThorLabs ECU Scanner

Configure a resonant scanner to link to a Thorlabs ECU1 device. This will configure the selected scanner to communicate with the ECU on the selected serial COM port.

Scanner Name: ResScanner

ECU COM Port: COM1

Select the resonant scanning system you defined in the [General ScanImage Settings](#) section. Only scanners that have a "Resonant" Scanner Type are available in the Scanner Name dropdown list.

Recap:

General ScanImage Settings

List all scanning systems to control on this microscope. For each, specify whether it is a resonant (resonant-galvo scan mirror pair) or linear (galvo-galvo scan mirror pair) system. A custom startup script and plugins can also be specified.

Scanning Systems

Scanner Type	Scanner Name	Delete
Resonant	ResScanner	X
Linear	LinScanner	X
		X

Select the ECU COM Port of USB to serial adapter that installs when USB cord was plugged in. Use the device manager to determine this.

 The Digital I/O device where triggers are wired is configured in the [Resonant Scanning system section](#).

Auxiliary Digital I/O DAQ

PXI1Slot3

Since the scanner control DAQ board is wired directly to the Thorlabs ECU, a separate board must be used. This can be found in NI MAX. (ex: 'PXI1Slot3').

Configure Thorlabs B-Scope 2

Select the ThorLabs BScope2 Setup section In the ScanImage Machine Configuration Editor.

ThorLabs BScope2 Setup

Configure a Thorlabs B-Scope 2. Select a resonant scanner to link to the Thorlabs ECU2 on the selected serial COM port. If USB mirror controllers are available select the appropriate COM port. If USB PMT's are available enter the VISA address or 6 digit serial numbers.

Scanner Name	<input type="text" value="ResScanner"/>	<input type="checkbox"/>
ECU COM Port	<input type="text" value="COM2"/>	<input type="checkbox"/>
Galvo-Galvo Flipper Mirror COM Port	<input type="text" value="COM3"/>	<input checked="" type="checkbox"/> Invert Mirror Position
Resonant-Galvo Flipper Mirror COM Port	<input type="text" value="COM1"/>	<input checked="" type="checkbox"/> Invert Mirror Position
Camera PMT Flipper Mirror COM Port	<input type="text" value="COM4"/>	<input checked="" type="checkbox"/> Invert Mirror Position

USB PMT Controller Serial Number or VISA Address	Delete
	X

Select the resonant scanning system you defined in the [General ScanImage Settings](#) section. Only scanners that have a "Resonant" Scanner Type are available in the Scanner Name dropdown list.

Recap:

General ScanImage Settings

List all scanning systems to control on this microscope. For each, specify whether it is a resonant (resonant-galvo scan mirror pair) or linear (galvo-galvo scan mirror pair) system. A custom startup script and plugins can also be specified.

Scanning Systems

Scanner Type	Scanner Name	Delete
Resonant	ResScanner	X
Linear	LinScanner	X
		X

Select the appropriate COM Ports for ECU, Galvo-Galvo Flipper mirror, Resonant-Galvo Flipper mirror, and Camera PMT Flipper mirror.

 Select the ECU COM Port of USB to serial adapter that installs when USB cord was plugged in. Use the device manager to determine this.

See [ThorLabs BScope2 Settings](#) for additional setup information.

 The FlexRIO Device is configured in the [Resonant Scanning system settings](#) section.

Signal Acquisition DAQ RIO0 (NI7961, NI5734) ▼

The available selections in the Signal Acquisition DAQ dropdown list are obtained from configured devices in NI MAX.

Configure Resonant Scanning System

Select the [Scanner Settings](#) section for the defined Resonant scanning system, In the ScanImage Machine Configuration Editor. The defined scanning system will appear in the section name.

Enter the additional Resonant and Galvo mirror configuration.

Resonant Mirror settings

Resonant Mirror

Zoom Control DAQ PXI1Slot3 ▼

Zoom Control Channel ID A00 ▼

Zoom Control Max Command (V) 4.95

Max Angular Range (optical deg pk-pk) 15

Zoom Amplitude Settling Time (sec) 0.5

Select the name of the DAQ device where Thorlabs ECU is plugged in (found in NI MAX. (ex: 'PXI1Slot3')) in the Zoom Control DAQ. Defaults are provided by the system. Change default entries as needed.

Galvo Mirror settings

Galvo Settings

Galvo Position Control DAQ PXI1Slot3 ▼

X Galvo (Optional)

Position Control AO Channel ID None ▼

Max Angular Range (optical deg pk-pk) 15

Command Scaling Factor (V/optical deg) 1

Park Angle (optical deg) -8

Y Galvo

Position Control AO Channel ID A01 ▼

Max Angular Range (optical deg pk-pk) 15

Command Scaling Factor (V/optical deg) 1

Park Angle (optical deg) -8

Defaults are provided by the system. Change default entries as needed.

